

SIIT-Chiba joint student workshop  
Presentation

# **Non-contact automatic measurement of hand joint range of motion**

Project Design : Prof. Ph.D Toshiya Nakaguchi

Group Members

Hironori Suzuki  
Haseeb Ali  
Yoshito Mikado  
Isura Dissanayake  
Pitchayagan Temniranrat

# Presentation Overview

- Measuring ROM
- Division of roles
- System overview
- Demo
- Evaluation
- Problems
- Future works

# Measuring ROM (Range of Motion)

Typical Measuring ROM of body joints require therapist to

- Adjust the pose of the target patient
- Use goniometer to measure the joint angles manually



# Measuring ROM's Problem

- Need physical contact
- Manual measurement is time consuming
- Large margin of error

Simple and accurate measurement of ROM using image analysis is highly required


→ Non-contact automatic measurement of hand joint ROM

# Division of Roles


- Program Flow Design
- GUI for image acquisition
- Hand Landmark Detection
- Angle calculation
- Image Annotations (Lines and Angle)
- Information organization and documentation

# System Overview

→ : save data  
→ : input data



# System Overview : Input image acquisition (GUI)



## CSV table for descriptions

[illegible]



# Joint coordinates detection

- MediaPipe (Our choice)

- Fan Zhang, Valentin Bazarevsky, Andrey Vakunov, Andrei Tkachenka, George Sung, Chuo-Ling Chang, & Matthias Grundmann. (2020). MediaPipe Hands: On-device Real-time Hand Tracking.
- Detect hand joints' position from single RGB camera

- Leap Motion

- “The Ultraleap’s LeapMotion Controller 2”, Ultraleap. Available: <https://www.ultraleap.com/product/leap-2/>. Accessed on : December 21, 2023.
- Used a specific device to detect hand joints' position
  - Includes of 2 cameras and an infrared sensor

# Hand LandMark Detection

## Image Input

Image Captured from GUI  
Stored into specified  
folder(to\_process)

## Mediapipe Hand Application

Fed Images to MediaPipe  
Library and obtained Hand  
Landmarks

## Output of Landmarks

All landmarks accumulated and  
saved as CSV, for further  
processing.  
(21 \* 3 rows and 36 \* 2 columns)

# Outputs After Landmark Detection

## Output CSV

	Joint_0_x	Joint_0_y	Joint_0_z	Joint_1_x	Joint_1_y	Joint_1_z
0	987.63393	352.86973	-0.00026	844.0468	366.84746	125.38770
1	1063.0320	331.87856	0.000623	1056.5705	387.81024	-147.061
2	1024.5969	816.79920	0.001062	934.58393	897.4641	-132.694
3	998.6937	754.44846	0.001456	969.66556	708.62743	-102.201
4	999.15206	494.07703	0.000776	908.89400	632.18846	4.9083875
5	1012.2698	508.61496	0.000299	976.96151	656.65002	8.3919677
6	1019.9899	476.59618	0.000711	928.2422	613.50368	8.9165228
7	1041.2881	471.32792	0.000776	904.13326	547.04459	26.482602
8	1033.4719	469.16406	0.000713	955.12144	612.19646	7.0469520
9	1031.5731	453.35639	0.001247	971.23329	593.37044	-1.94412
10	1032.2495	477.24573	0.000714	950.79906	618.32219	11.359963
11	1033.8219	469.31382	0.000732	947.85575	614.49696	0.3320946
12	1030.3405	394.6049	0.00053	906.74016	420.92581	-17.5075
13	1009.0090	383.1775	0.001807	947.45458	435.61542	-155.352

65\_processed



71\_processed



63\_processed



# Image Annotations

- Used OpenCV - cv2 Library for annotating images
- Outputs from Angle Calculation Program were read in
- Calculation of intersection and projection of points required for some lines.


## Angle.csv

Pose	Joint name	Angle
1	R -Wrist -	62.85358984532308
3	R -Wrist -	67.72033
5	R -Wrist -	20.88606076142549
7	R -Wrist -	29.520442921685333
9	R -Thumb -	77.55556733035968
11	R -Thumb -	8.072896312053095
13	R -Thumb -	57.929086172877135
15	R -Thumb -	13.785923884258336
17	R -Thumb -	46.21461027562268

Read In

## Lines Calculation

running_n_index	capture_p	sample_p	out_pic_n	basic_axis	basic_axis	moving_a	moving_a	line1_o	line1_o_c	line1_o_s	line1_v	line1_v_c	line1_v_s	line2_o
0	2	1	1	1		0	9	0	0	0	9	0	0	0
1	4	3	3	3		0	9	0	2	2	9	2	2	0
2	6	5	5	5		0	9	0	4	4	9	4	4	0
3	8	7	7	7		0	9	0	6	6	9	6	6	0
4	10	8	9	9	5	8	1	4 cutting_pc	8	9	8	8	9 cutting_pc	
5	12	9	11	11	5	8	1	4 proj(0,[5,8	9	11	8	9	11 proj(0,[1,4	
6	14	10	13	13	5	8	1	4 cutting_pc	10	13	8	10	13 cutting_pc	
7	16	11	15	15	5	8	1	4 proj(0,[5,8	11	15	8	11	15 proj(0,[1,4	






# Output Examples

## 1) CSV Table List

Pose	Joint name	Angle
29_0	R-Finger_0 - Flexion PIP	49.6985
29_1	R-Finger_1 - Flexion PIP	53.0159
29_2	R-Finger_2 - Flexion PIP	53.6065
29_3	R-Finger_3 - Flexion PIP	46.3057
31_0	R-Finger_0 - Extension PIP	10.1214
31_1	R-Finger_1 - Extension PIP	10.7499
31_2	R-Finger_2 - Extension PIP	10.2342
31_3	R-Finger_3 - Extension PIP	6.93028
33_0	R-Finger_0 - Flexion DIP	45.9127
33_1	R-Finger_1 - Flexion DIP	46.2441
33_2	R-Finger_2 - Flexion DIP	54.2222
33_3	R-Finger_3 - Flexion DIP	62.3758
35_0	R-Finger_0 - Extension DIP	9.18794
35_1	R-Finger_1 - Extension DIP	7.61211
35_2	R-Finger_2 - Extension DIP	11.0205
35_3	R-Finger_3 - Extension DIP	8.78998
37	L-Wrist - Flexion	61.3588
39	L-Wrist - Extension	71.4339
41	L-Wrist - Radial Deviation	20.6558
43	L-Wrist - Ulnar Deviation	29.3692
45	L-Thumb - Radial Abduction	71.4531
47	L-Thumb - Ulnar Abduction	9.88693
49	L-Thumb - Palmer Abduction	57.1438
51	L-Thumb - Palmer Abduction	13.2684
53	L-Thumb - Flexion MCP	50.4373
55	L-Thumb - Extension MCP	2.90352
57	L-Thumb - Flexion IP	65.6583

## 2) Measured Angles



# Demo







# Evaluation

Our angles measured manually may be not accurate because we do not know how to measure rightly

Mean squared error value is **14.32 degrees** (middle-ring-little fingers not included)



# Problems and Future Works

- Joint landmarks detection
  - Media Pipe : Inaccurate joint landmark detection in case of occlusion
  - Leap Motion : More accurate results, but difficult coding on python (other parts of the system are on python)
- Axis - Corresponding Landmarks Decision
  - Many cases to consider without common factor → have to be manually prepared
- Future works :
  - Using multi-camera to improve the performance of joint detection
  - Adding a retake function if the user are not satisfied with the captured image